

PRODUCT RANGE SPECIFICATION

H65240 – Instrument Trolley with Bars

DOCUMENT SCOPE Product range	DOCUMENT REFERENCE H65240{PSS}	DATA STATUS Published catalogue data	GENERATED ON 14 Sep 2026
---------------------------------	-----------------------------------	---	-----------------------------

01 Product description and specification

Product description	Instrument trolley with retaining bars for organising and moving instruments in healthcare treatment and procedure areas. The trolley is designed to ensure that instruments are easy to reach during operating procedures ensuring efficiency within the theater. Manufactured from Grade 304 stainless steel, the instrument trolley is easy to clean and maintain. Brushed polish stainless steel finish is aesthetic but also adds to the ease of cleaning of the trolley. Our H65240 Instrument Trolley is built to last, featuring a strong substructure made from durable stainless steel tubular frame. The theater instrument trolley is available with either two or three shelves or trays. The shelves are reinforced with underbracing which prevents flex in the shelf. The underbracing support also ensures that the shelves are rigid and strong. The trolley includes 4 × 125 mm castors which make the instrument trolley incredibly easy to move and position within the operating theater. The castors are also lockable which allows the instrument trolley to be fixed in place. This prevents the trolley from being accidentally bumped and moved during procedures. Construction and features • Surgical instrument trolley. • The instrument trolley is designed to be easy to position and ergonomic to use. • Manufactured from Grade 304 stainless steel ensuring high sanitary standards are exceeded. • The substructure is manufactured from a durable stainless steel tubular frame. • Available with either two or three shelves. • Reinforced underbracing ensures that the shelves are rigid and do not flex. • All stainless steel surfaces have an aesthetic brushed polish stainless steel finish. • Features 4 × 125 mm lockable castors as standard.
---------------------	--

Short Form Specification	SPECIFICATION USE Use this description when specifying this product range. Basterfield H65240. Instrument Trolley with Bars. Length: 490 mm, 600 mm, 900 mm, 1200 mm. Shelf Tiers: 2, 3. The trolley includes 4 × 125 mm castors which make the instrument trolley incredibly easy to move and position within the operating theater. Features 4 × 125 mm lockable castors as standard.
--------------------------	---

02 Models and configurations

MODEL	DESCRIPTION	PUBLISHED CONFIGURATION
H6524004902T	H6524004902T Two-Tier Instrument Trolley 490mm Length Shelf Tiers: 2. Length: 490.	Length: 490 mm Shelf Tiers: 2

MODEL	DESCRIPTION	PUBLISHED CONFIGURATION
H6524004903T	H6524004903T Three-Tier Instrument Trolley 490mm Length	Length: 490 mm Shelf Tiers: 3
	Shelf Tiers: 3. Length: 490.	
H6524006002T	H6524006002T Two-Tier Instrument Trolley 600mm Length	Length: 600 mm Shelf Tiers: 2
	Shelf Tiers: 2. Length: 600.	
H6524006003T	H6524006003T Three-Tier Instrument Trolley 600mm Length	Length: 600 mm Shelf Tiers: 3
	Shelf Tiers: 3. Length: 600.	
H6524009002T	H6524009002T Two-Tier Instrument Trolley 900mm Length	Length: 900 mm Shelf Tiers: 2
	Shelf Tiers: 2. Length: 900.	
H6524009003T	H6524009003T Three-Tier Instrument Trolley 900mm Length	Length: 900 mm Shelf Tiers: 3
	Shelf Tiers: 3. Length: 900.	
H6524012002T	H6524012002T Two-Tier Instrument Trolley 1200mm Length	Length: 1200 mm Shelf Tiers: 2
	Shelf Tiers: 2. Length: 1200.	
H6524012003T	H6524012003T Three-Tier Instrument Trolley 1200mm Length	Length: 1200 mm Shelf Tiers: 3
	Shelf Tiers: 3. Length: 1200.	


03 Application information


Intended applications	Healthcare Utility Rooms
Industries	Hospitals & Clinics

04 Supply and support


Typical lead time	10–15 working days
Warranty	12 months

Product drawings and images


**PRODUCT IMAGE**



PRODUCT IMAGE



PRODUCT IMAGE




PRODUCT IMAGE




PRODUCT IMAGE



PRODUCT IMAGE



PRODUCT IMAGE



PRODUCT IMAGE

VERIFICATION REQUIRED Published values are provided for initial specification. Confirm the selected model, critical dimensions, site services, clearances, and the latest document revision with Basterfield before procurement, manufacture, or installation. E&OE. Specifications are subject to change without notice.