

PRODUCT RANGE SPECIFICATION

B36658 – Stainless Steel Knee Activated Wash Hand Basin

DOCUMENT SCOPE Product range	DOCUMENT REFERENCE B36658{PSS}	DATA STATUS Published catalogue data	GENERATED ON 14 Sep 2026
---------------------------------	-----------------------------------	---	-----------------------------

01 Product description and specification

Product description	The B36658 knee-operated wash hand basin has an Easy Push front panel for operation with either knee. A gooseneck spout provides aerated water, and the stainless steel basin is supported by an easy-mount bracket with multiple fixing holes. Choose a demand valve for flow while the panel is pressed, or a metered valve for a preset flow period after activation. An extended stainless steel splashback is optional. Confirm the mounting position, valve option and access requirements for the installation.
---------------------	---

Short Form Specification	SPECIFICATION USE Use this description when specifying this product range. Basterfield B36658. Stainless steel knee-operated wash basin; Easy Push front panel; aerated gooseneck spout; wall mounting bracket. Demand or metered valve as selected; extended splashback optional.
--------------------------	--

02 Models and configurations

MODEL	DESCRIPTION	PUBLISHED CONFIGURATION
B36658(D)	KNEE OPERATED WASH BASIN WITH DEMAND VALVE ACTIVATION	Demand / Metered Tap: Demand Tap Extended Splashback: No
B36658(DE)	KNEE OPERATED WASH BASIN WITH EXTENDED SPLASH BACK AND DEMAND VALVE ACTIVATION	Demand / Metered Tap: Demand Tap Extended Splashback: Yes
B36658(ME)	KNEE OPERATED WASH BASIN WITH EXTENDED SPLASH BACK AND METERED VALVE ACTIVATION	Demand / Metered Tap: Metered Tap Extended Splashback: Yes
B36658M	B36658M Knee Activated Wash Basin fitted with metered tap Extended Splashback: No. Demand / Metered Tap: Metered Tap.	Demand / Metered Tap: Metered Tap Extended Splashback: No

03 Application information

Intended applications	Hands-Free Handwashing (Knee/Elbow/Foot Operated)
Industries	Government & Public Sector, Hospitals & Clinics, Research Laboratories

04 Supply and support

Typical lead time	10–15 working days
Warranty	12 months

Product drawings and images**PRODUCT IMAGE**

VERIFICATION REQUIRED Published values are provided for initial specification. Confirm the selected model, critical dimensions, site services, clearances, and the latest document revision with Basterfield before procurement, manufacture, or installation. E&OE. Specifications are subject to change without notice.